

Реагенты для клинического анализа

Технические характеристики

По вопросам продаж и поддержки обращайтесь:

Алматы (727)345-47-04
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Саранск (8342)22-96-24
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +(727)345-47-04

Беларусь +(375)257-127-884

Узбекистан +998(71)205-18-59

Киргизия +996(312)96-26-47

эл.почта: cmf@nt-rt.ru || сайт: <https://chem-lab.nt-rt.ru/>

Azur-Eosine-Methylenblue solution according Giemsa
8 g Giemsa's stain / l glycerol/methanol (1/1)
Microbiologic colourations

CL04.0502

Physical

- **Form:** liquid
- **Colour:** blue
- **Melting point:** -98°C
- **Boiling point:** 65°C
- **Flash point:** 11°C

- **Density:** 1,02 g/ml
- **Storage temp:** RT

Art. Nr.	Pack	Mod	Stock
CL04.0502.1000	1 l	GVB	0

CL04.0201

BK-Auramine solution (A)

1 g C₁₇H₂₁N₃.HCl + 25.5 g C₆H₅OH / l H₂O

Art. Nr.	Pack	Mod	Stock
CL04.0201.1000	1 l	PE	0

CL04.0202

BK-Decolouration solution (B)

5 g NaCl + 300 ml H₂O + 5 ml HCl 37% + 700 ml C₂H₅OH

Art. Nr.	Pack	Mod	Stock
CL04.0202.1000	1 l	PE	0

CL04.0203

BK-Methylene blue solution (C)

20 g C₁₆H₁₈CIN₃S.3H₂O / l H₂O

Casnr.: 61-73-4

Art. Nr.	Pack	Mod	Stock
CL04.0203.1000	1 l	PE	0

CL04.0204
Bouin solution
 100 ml HCOH 37% + 6 g C6H3N3O7 + 50 ml CH3COOH / H2O

Art. Nr.	Pack	Mod	Stock
CL04.0204.1000	1 l	PEB	2

CL04.0301
Cresyl violet acetate
 C18H15N3O3
 Casnr.: 10510-54-0

Art. Nr.	Pack	Mod	Stock
CL04.0301.0005	5 g	GVB	0

CL04.0703
Crystal Violet
 C25H30ClN3
 Casnr.: 548-62-9

Art. Nr.	Pack	Mod	Stock	
CL04.0703.0025	25 g	GVB	0	
CL04.0703.0100	100 g	GVB	1	
CL04.0703.9025	25 kg	CBOT	0	Custom

CL04.1102
Crystal Violet , 1% w/v in Acetic acid 99-100%
 10 g C25H30ClN3 / l CH3COOH

Art. Nr.	Pack	Mod	Stock
CL04.1102.0100	100 ml	GVB	8

CL04.0504
Eosine Blue (Eosine B)
C₂₀H₆Br₂N₂Na₂O₉
Casn.: 548-24-3

Art. Nr.	Pack	Mod	Stock
CL04.0504.0025	25 g	GVB	3

CL04.0501
Eosine-Methylene blue according May-Grünwald
2.5 g eosine-methyleenblauw / l CH₃OH

Art. Nr.	Pack	Mod	Stock
CL04.0501.0250	250 ml	GVB	3
CL04.0501.1000	1 l	GVB	5

CL04.0506
Eosine Yellow (Eosine Y)
C₂₀H₆Br₄Na₂O₅
Casn.: 17372-87-1

Art. Nr.	Pack	Mod	Stock
CL04.0506.0025	25 g	GVB	19

CL04.0505
Erythrosine 0.5 g/l in Propanol-2 (iso-Propanol) a.r.
0.5 g C₂₀H₆I₄Na₂O₅ / l C₃H₈O

Art. Nr.	Pack	Mod	Stock	
CL04.0505.5000	5 l	PE	15	Custom

CL04.0503
Eukitt
n.a.

Art. Nr.	Pack	Mod	Stock	
CL04.0503.0100	100 ml	GVB	0	Custom
CL04.0503.0500	500 ml	GVB	1	Custom

CL04.0602

Fuchsine 2% Ethanol solution

20 g C₂₀H₂₀CIN₃ / l C₂H₅OH

Art. Nr.	Pack	Mod	Stock	
CL04.0602.0250	250 ml	GVB	2	
CL04.0603				
Fuchsine acc. Ziehl solution				
10 g C ₂ OH ₂ OCIN ₃ + 50 g C ₆ H ₅ OH + 100 ml C ₂ H ₅ OH / l H ₂ O				

Art. Nr.	Pack	Mod	Stock	
CL04.0603.0100	100 ml	GVB	8	
CL04.0603.1000	1 l	GVB	19	
CL04.0603.2500	2.5 l	GVB	3	

CL04.0801

Hemalein acc. Mayer solution

20 g KAl(SO₄)₂.12H₂O + 2 g C₁₆H₁₄O₆ + 20 ml CH₃COOH + 90 ml CH₃OH / l H₂O

Art. Nr.	Pack	Mod	Stock	
CL04.0801.1000	1 l	GVB	0	
CL04.0801.2500	2.5 l	GVB	0	

CL04.0901

Immersion-oil

25 à 50% C₁₄H₁₂O₂ (Benzyl benzoate)

Casnr.: 120-51-4

Art. Nr.	Pack	Mod	Stock
CL04.0901.0100	100 ml	GVB	16
CL04.0901.0250	250 ml	GVB	4

CL04.1202

Lactophenol acc. Amann solution

200 g C₆H₅OH + 200 g C₃H₆O₃ + 400 ml C₃H₈O₃ / l H₂O

Art. Nr.	Pack	Mod	Stock
CL04.1202.0125	125 ml	GVB	3

CL04.1204

Lactophenol - Aniline blue solution

200 g C₆H₅OH + 200 g C₃H₆O₃ + 400 ml C₃H₈O₃ + 0.5 g C₃₂H₂₅N₃Na₂O₉S₃ / l H₂O

Art. Nr.	Pack	Mod	Stock
CL04.1204.0100	100 ml	GVB	1

CL04.1206

Lugol solution

3.4 g I₂ + 6.3 g KI / l H₂O

Art. Nr.	Pack	Mod	Stock	
CL04.1206.1000	1 l	GVB	4	Custom

Lugol solution

50 g I₂ + 100 g KI / l H₂O

Microbiologic colourations

CL04.1201

Physical

- **Form:** liquid
- **Colour:** brown
- **Melting point:** 0°C
- **Boiling point:** 100°C
- **Flash point:** -

- **Density:** 1,20 g/ml
- **Storage temp:** RT

Art. Nr.	Pack	Mod	Stock
CL04.1201.0500	500 ml	GVB	12
CL04.1201.1000	1 l	GVB	0

CL04.1203

Lugol - PVP solution

5 g I₂ + 10 g KI + 88 g (C₆H₉NO)_nxl / l H₂O

Art. Nr.	Pack	Mod	Stock	
CL04.1203.0100	100 ml	GVB	0	Custom

CL04.1302

Methylene blue alkalic solution acc. Loeffler

16.2 ml KOH 0.1 mol/l + 91 ml C₂H₅OH + 3.65 g C₁₆H₁₈CIN₃S / l H₂O

Art. Nr.	Pack	Mod	Stock	
CL04.1302.0100	100 ml	PE	1	Custom
CL04.1302.0250	250 ml	PE	0	

CL04.1301

Phenol - Methylene blue solution

10 g C₆H₅OH + 100 g C₂H₅OH + 15 g C₁₆H₁₈CIN₃S / l H₂O

Art. Nr.	Pack	Mod	Stock
CL04.1301.0125	125 ml	GVB	5
CL04.1301.1000	1 l	GVB	0
CL04.1301.2500	2.5 l	GVB	0

CL04.0702

Phenol - Crystal Violet solution

20 g C₆H₅OH + 100 ml C₂H₅OH + 10 g C₂₅H₃₀CIN₃ / l H₂O

Art. Nr.	Pack	Mod	Stock	
CL04.0702.1000	1 l	GVB	0	Custom
CL04.0702.2500	2.5 l	GVB	0	

CL04.0704

Phenol - Crystal Violet solution

4 g C₆H₅OH + 10 g C₂₅H₃₀CIN₃ / l H₂O

Art. Nr.	Pack	Mod	Stock
CL04.0704.1000	1 l	GVB	0

CL04.2001

Phenol - Toluidine blue O solution

30 g C₆H₅OH + 360 ml CH₃OH + 100 ml C₂H₅OH + 5.9 g C₁₅H₁₆CIN₃S / l H₂O

Art. Nr.	Pack	Mod	Stock
CL04.2001.0100	100 ml	GVB	0

CL04.0604

Phenolphthalein in Ethanol

0.5 g C₂₀H₁₄O₄ / l C₂H₅OH

Art. Nr.	Pack	Mod	Stock
CL04.0604.0100	100 ml	GVB	0

CL04.1601

Ponceau S red solution

2 g Ponceau S red / l 3% C₂HCl₃O₂ in H₂O

Art. Nr.	Pack	Mod	Stock
CL04.1601.0250	250 ml	PE	1

CL04.1901

Safranine O acc. Gram-Hucker solution

4 g C₂₀H₁₉CIN₄ / l C₂H₅OH 20%

Art. Nr.	Pack	Mod	Stock
CL04.1901.0250	250 ml	GVB	4
CL04.1901.1000	1 l	GVB	0

CL04.1902

Safranine O solution

1.8 g C₂₀H₁₉CIN₄ / l H₂O

Casnr.: 477-73-6

Art. Nr.	Pack	Mod	Stock
CL04.1902.1000	1 l	GVB	3

По вопросам продаж и поддержки обращайтесь:

Алматы (727)345-47-04
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Саранск (8342)22-96-24
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +(727)345-47-04

Беларусь +(375)257-127-884

Узбекистан +998(71)205-18-59

Киргизия +996(312)96-26-47

эл.почта: cmf@nt-rt.ru || сайт: <https://chem-lab.nt-rt.ru/>