

Стандарты электродов для измерений

Технические характеристики

По вопросам продаж и поддержки обращайтесь:

Алматы (727)345-47-04
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Саранск (8342)22-96-24
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +(727)345-47-04

Беларусь +(375)257-127-884

Узбекистан +998(71)205-18-59

Киргизия +996(312)96-26-47

эл.почта: cmf@nt-rt.ru || сайт: <https://chem-lab.nt-rt.ru/>

Hydrogen chloride solution 500000 µS/cm (Conductivity standard)
LF Standard 500000 µS/cm @ /25°C (± 1%)
For laboratory use, LF measurements
Traceability: Reference standards acc. NIST SRM 999b

CL03.1149

CAS Nr:7647-01-0

Physical

- **Form:** liquid
- **Colour:** colourless
- **Melting point:** 0°C
- **Boiling point:** 100°C
- **Flash point:** -

- **Density:** 1,00 g/ml
- **Mol Weight:** 36.46 g/mol
- **Storage temp:** 2°C to 25°C

Art. Nr.	Pack	Mod	Stoc k
CL03.1149.0500	500 ml	PE	0

CL03.1148

Hydrogen chloride solution 450000 µS/cm (Conductivity standard)
LF Standard 450000 µS/cm @ /25°C (± 1%)
Casnr.: 7647-01-0

Art. Nr.	Pack	Mod	Stoc k
CL03.1148.0500	500 ml	PE	0

CL03.1147

Hydrogen chloride solution 350000 µS/cm (Conductivity standard)
LF Standard 350000 µS/cm @ /25°C (± 1%)
Casnr.: 7647-01-0

Art. Nr.	Pack	Mod	Stoc k
CL03.1147.0500	500 ml	PE	0

CL03.1146
Potassium chloride solution 300000 µS/cm (Conductivity standard)
 LF Standard 300000 µS/cm @ 25°C (± 1%)
 Casnr.: 7447-40-7

Art. Nr.	Pack	Mod	Stoc k
CL03.1146.0500	500 ml	PE	0

CL03.1145
Potassium chloride solution 200000 µS/cm (Conductivity standard)
 LF Standard 200000 µS/cm @ 25°C (± 1%)
 Casnr.: 7447-40-7

Art. Nr.	Pack	Mod	Stoc k
CL03.1145.0500	500 ml	PE	0

CL03.1145.9010	10 l	PE	0	Custom
----------------	------	----	---	--------

CL03.1144
Potassium chloride solution 150000 µS/cm (Conductivity standard)
 LF Standard 150000 µS/cm @ 25°C (± 1%)
 Casnr.: 7447-40-7

Art. Nr.	Pack	Mod	Stoc k
CL03.1144.0500	500 ml	PE	8

CL03.1104
Potassium chloride 1 mol/l (Conductivity standard)
 LF Standard 102.1 mS/cm @ 20°C - 111.8 mS/cm @ 25°C (± 1%)
 Casnr.: 7447-40-7

Art. Nr.	Pack	Mod	Stoc k
CL03.1104.1000	1 l	PE	46
CL03.1104.9010	10 l	PE	0

CL03.1143
Potassium chloride solution 100000 µS/cm (Conductivity standard)
 LF Standard 100000 µS/cm @ 25°C (± 1%)
 Casnr.: 7447-40-7

Art. Nr.	Pack	Mod	Stoc k
CL03.1143.0500	500 ml	PE	4

CL03.1142
Potassium chloride solution 50000 µS/cm (Conductivity standard)
 LF Standard 50000 µS/cm @ 25°C (± 1%)
 Casnr.: 7447-40-7

Art. Nr.	Pack	Mod	Stoc k
CL03.1142.0500	500 ml	PE	0

CL03.1126
Potassium chloride solution 46 mS/cm (Conductivity standard)
 LF Standard 46 mS/cm @ 20°C - 50 mS/cm @ 25°C (± 1%)
 Casnr.: 7447-40-7

Art. Nr.	Pack	Mod	Stoc k
CL03.1126.0250	250 ml	PE	0

CL03.1126.0500	500 ml	PE	1
----------------	-----------	----	---

CL03.1124
Potassium chloride 0.2 mol/l (Conductivity standard)
 LF Standard 22.4 mS/cm @ 20°C - 24.7 mS/cm @ 25°C (± 1%)
 Casnr.: 7447-40-7

Art. Nr.	Pack	Mod	Stoc k
CL03.1124.0500	500 ml	PE	1

CL03.1124.1000	1 l	PE	8
----------------	-----	----	---

CL03.1141
Potassium chloride solution 20000 µS/cm (Conductivity standard)
 LF Standard 20000 µS/cm @ 25°C (± 1%)
 Casnr.: 7447-40-7

Art. Nr.	Pack	Mod	Stoc k
CL03.1141.0500	500 ml	PE	1

CL03.1108

Potassium chloride 0.1 mol/l (Conductivity standard)

LF Standard 11.67 mS/cm @ 20°C - 12.88 mS/cm @ 25°C (± 1%)

Casnr.: 7447-40-7

Art. Nr.	Pack	Mod	Stoc k
CL03.1108.0500	500 ml	PE	40

CL03.1108.1000	1 l	PE	47
----------------	-----	----	----

CL03.1108.9010	10 l	PE	11
----------------	------	----	----

CL03.1140

Potassium chloride solution 10000 µS/cm (Conductivity standard)

LF Standard 10000 µS/cm @ 25°C (± 1%)

Casnr.: 7447-40-7

Art. Nr.	Pack	Mod	Stoc k
CL03.1140.0500	500 ml	PE	3

CL03.1139

Conductivity standard 5000 µS/cm (Potassium chloride solution)

LF Standard 5000 µS/cm @ 25°C (± 1%)

Casnr.: 7447-40-7

Art. Nr.	Pack	Mod	Stoc k
CL03.1139.0500	500 ml	PE	2

CL03.1103

Potassium chloride 0.03 mol/l (Conductivity standard)

LF Standard 3660 $\mu\text{S/cm}$ @ 20°C - 4098 $\mu\text{S/cm}$ @ 25°C ($\pm 1\%$)

Casnr.: 7447-40-7

Art. Nr.	Pack	Mod	Stoc k
CL03.1103.1000	1 l	PE	0

CL03.1103.9010	10 l	PE	0
----------------	------	----	---

CL03.1102

Potassium chloride 0.02 mol/l (Conductivity standard)

LF Standard 2542 $\mu\text{S/cm}$ @ 20°C - 2765 $\mu\text{S/cm}$ @ 25°C ($\pm 1\%$)

Casnr.: 7447-40-7

Art. Nr.	Pack	Mod	Stoc k
CL03.1102.1000	1 l	PE	2

CL03.1102.9010	10 l	PE	0
----------------	------	----	---

CL03.1101

Potassium chloride 0.01 mol/l (Conductivity standard)

LF Standard 1278 $\mu\text{S/cm}$ @ 20°C - 1413 $\mu\text{S/cm}$ @ /25°C ($\pm 1\%$)

Casnr.: 7447-40-7

Art. Nr.	Pack	Mod	Stoc k
CL03.1101.0500	500 ml	PE	187
CL03.1101.1000	1 l	PE	95
CL03.1101.9010	10 l	PE	0

CL03.1120

Potassium chloride solution 1100 µS/cm (Conductivity standard)

LF Standard 1100 µS/cm @ 25°C (± 1%)

Casnr.: 7447-40-7

Art. Nr.	Pack	Mod	Stoc k
CL03.1120.0500	500 ml	PE	2

CL03.1138

Potassium chloride solution 1000 µS/cm (Conductivity standard)

LF Standard 1000 µS/cm @ 25°C (± 1%)

Casnr.: 7447-40-7

Art. Nr.	Pack	Mod	Stoc k
CL03.1138.0500	500 ml	PE	3

Art. Nr.	Pack	Mod	Stoc k
CL03.1121.0500	500 ml	PE	48

CL03.1153

Potassium chloride solution 445 µS/cm (Conductivity standard)

LF Standard 445 µS/cm @ 25°C (± 1%)

Casnr.: 7447-40-7

Art. Nr.	Pack	Mod	Stoc k
CL03.1153.1000	1 l	PE	0

CL03.1137

Potassium chloride solution 200 µS/cm (Conductivity standard)

LF Standard 200 µS/cm @ 25°C (± 1%)

Casnr.: 7447-40-7

Art. Nr.	Pack	Mod	Stoc k
CL03.1137.0500	500 ml	PE	1

CL03.1133

Potassium chloride 0.001 mol/l (Conductivity standard)

LF Standard 133 µS/cm @ 20°C - 147 µS/cm @ 25°C (± 1%)

Casnr.: 7447-40-7

Art. Nr.	Pack	Mod	Stoc k
CL03.1133.0500	500 ml	PE	37

CL03.1119

Potassium chloride solution 100 µS/cm (Conductivity standard)

LF Standard 100 µS/cm @ 25°C (± 1%)

Casnr.: 7447-40-7

Art. Nr.	Pack	Mod	Stoc k
CL03.1119.0500	500 ml	PE	8

CL03.1150

Potassium chloride solution 84 µS/cm (Conductivity standard)

LF Standard 84 µS/cm @ 25°C (± 1%)

Casnr.: 7447-40-7

Art. Nr.	Pack	Mod	Stoc k
CL03.1150.0500	500 ml	PE	49

CL03.1134

Potassium chloride solution 50 µS/cm (Low Cond. standard)

LF Standard 50 µS/cm @ 25°C (± 1%) (Keep Cool !)

Casnr.: 7447-40-7

Art. Nr.	Pack	Mod	Stoc k
CL03.1134.0500	500 ml	PE	1

CL03.1128

Potassium chloride solution 25 μ S/cm (Low Cond. standard)

LF Standard 25 μ S/cm @ 25°C ($\pm$ 1%) (Keep Cool !)

Casnr.: 7447-40-7

Art. Nr.	Pack	Mod	Stoc k
CL03.1128.0500	500 ml	PE	2

CL03.1128.1000	1 l	PE	0
----------------	-----	----	---

CL03.1132

Potassium chloride solution 25 μ S/cm (Low Cond. standard)

LF Standard 25 μ S/cm @ 20°C ($\pm$ 1%) (Keep Cool !)

Casnr.: 7447-40-7

Art. Nr.	Pack	Mod	Stoc k
CL03.1132.1000	1 l	PE	0

CL03.1136

Potassium chloride solution 20 μ S/cm (Ultralow Cond. standard)

LF Standard 20 μ S/cm @ 25°C ($\pm$ 1%) (Keep Cool !)

Casnr.: 7447-40-7

Art. Nr.	Pack	Mod	Stoc k
CL03.1136.0500	500 ml	PE	0

CL03.1135

Potassium chloride solution 10 μ S/cm (Ultralow Cond. standard)

LF Standard 10 μ S/cm @ 25°C ($\pm$ 1%) (Keep Cool !)

Casnr.: 7447-40-7

Art. Nr.	Pack	Mod	Stoc k
CL03.1135.0500	500 ml	PE	3

CL03.1127

Potassium chloride solution 5 μ S/cm (Ultralow Cond. standard)

LF Standard 5 μ S/cm /@ 25°C ($\pm$ 1%) (Keep Cool !)

Casnr.: 7447-40-7

Art. Nr.	Pack	Mod	Stoc k
CL03.1127.0500	500 ml	PE	6

Art. Nr.	Pack	Mod	Stoc k
CL03.1805.0500	500 ml	PE	2

CL03.1806

Redox standard solution +470 mV

(470 $\pm$ 5) mV - 20°C

Art. Nr.	Pack	Mod	Stoc k	
CL03.1806.9010	10 l	PE	0	Custom

CL03.1801

Redox standard solution +430 mV
(430 ± 5) mV 25°C

Art. Nr.	Pack	Mod	Stoc k	
CL03.1801.0500	500 ml	PE	1	
CL03.1801.1000	1 l	PE	0	

CL03.1804

Redox standard solution +358 mV
(358 ± 5) mV - 25°C

Art. Nr.	Pack	Mod	Stoc k	
CL03.1804.0500	500 ml	PE	5	
CL03.1804.1000	1 l	PE	0	

CL03.1802

Redox standard solution +183 mV
(183 ± 5) mV - 25°C

Art. Nr.	Pack	Mod	Stoc k
CL03.1802.0500	500 ml	PE	0
CL03.1802.1000	1 l	PE	0

CL03.1803

Redox standard solution +124 mV
(124 ± 5) mV - 25°C

Art. Nr.	Pack	Mod	Stoc k	
CL03.1803.0500	500 ml	PE	2	
CL03.1803.1000	1 l	PE	0	
CL03.1803.9010	10 l	PE	0	Custom

Art. Nr.	Pack	Mod	Stoc k	
CL03.0122.1000	1 l	PE/OD	22	
CL03.0122.9010	10 l	PE/OD	0	Custom

CL03.0905

Ammonium sulfate 4 mol/l - ISE buffer
528 g (NH₄)₂SO₄ / l H₂O
Casn.: 7783-20-2

Art. Nr.	Pack	Mod	Stoc k
CL03.0905.0500	500 ml	PE	0

CL03.0912
Ammonium sulfate 2 mol/l electrode filling solution - ISE buffer
 264 g (NH₄)₂SO₄ / l H₂O
 Casnr.: 7783-20-2

Art. Nr.	Pack	Mod	Stoc k
CL03.0912.0100	100 ml	PE	0

CL02.1159
Copper(II)-di-Ammonium EDTA 0.1 mol/l - ISE buffer
 38.784 g Cu(NH₄)₂C₁₀H₁₂N₂O₈ / l H₂O
 Casnr.: 67989-88-2

Art. Nr.	Pack	Mod	Stoc k
CL02.1159.0500	500 ml	PE	1

CL02.1197
Copper(II)-di-Ammonium EGTA 0.05 mol/l - ISE buffer
 27.21 g CuK₂(NH₄)₂C₁₄H₂₀N₂O₁₀ / l H₂O

Art. Nr.	Pack	Mod	Stoc k
CL02.1197.1000	1 l	PE	0

CL03.0908

Lithium chloride 4 mol/l - ISE buffer

170 g LiCl / l H₂O

Casnr.: 7447-41-8

Art. Nr.	Pack	Mod	Stoc k
CL03.0908.0500	500 ml	PE	0

CL01.1444

Nitrate standard solution kit for ISE measurements

0 - 1 - 10 - 100 & 1000 mg/l NO₃ in 1% KAl(SO₄)₂·12H₂O sol.

Art. Nr.	Pack	Mod	Stoc k
CL01.1444.0001	1 pc	PE	0

CL03.1123

Nitrate ISE interference buffer solution - ISE buffer

17.32g Al₂(SO₄)₃·18H₂O+3.43g Ag₂SO₄+1.28g H₃BO₃+2.52g H₃NO₃S / l H₂O

Art. Nr.	Pack	Mod	Stoc k
CL03.1123.0500	500 ml	GVB	0

CL03.0902

Potassium chloride 4 mol/l - ISE buffer

305 g KCl / l H₂O

Casnr.: 7447-40-7

Art. Nr.	Pack	Mod	Stoc k
CL03.0902.0500	500 ml	PE	0

CL03.1117
Potassium chloride 3.5 mol/l+Aluminium(III) sulfate 50 g/l solution - ISA buffer - ISE buffer
 250 g KCl + 50 g AlK(SO4)2.18H2O / l H2O

Art. Nr.	Pack	Mod	Stoc k
CL03.1117.0500	500 ml	PE	0

CL03.1117.1000	1 l	PE	0
----------------	-----	----	---

CL03.0906
Potassium nitrate 2.5 mol/l - ISE buffer
 253 g KNO3 / l H2O
 Casnr.: 7757-79-1

Art. Nr.	Pack	Mod	Stoc k
CL03.0906.0500	500 ml	PE	0

CL03.1407
Sodium acetate buffer pH 4.7 - ISE buffer
 132 g CH3COONa.3H2O + 86 ml CH3COOH / l H2O
 Casnr.: 6131-90-4

Art. Nr.	Pack	Mod	Stoc k
CL03.1407.1000	1 l	PE	0

CL03.0904
Sodium chloride 5 mol/l - ISE buffer
 292 g NaCl / l H₂O
 Casnr.: 7647-14-5

Art. Nr.	Pack	Mod	Stoc k
CL03.0904.0500	500 ml	PE	0

CL03.0901
Sodium hydroxide 10 mol/l - ISE buffer
 400 g NaOH / l H₂O
 Casnr.: 1310-73-2

Art. Nr.	Pack	Mod	Stoc k
CL03.0901.0500	500 ml	PE	0

CL03.0903
Sodium nitrate 5 mol/l - ISE buffer
 425 g NaNO₃ / l H₂O
 Casnr.: 7631-99-4

Art. Nr.	Pack	Mod	Stoc k
CL03.0903.0500	500 ml	PE	1

CL03.0907

Sulfuric acid 2 mol/l - ISE buffer

196 g H₂SO₄ / l H₂O

Casnr.: 7664-93-9

Art. Nr.	Pack	Mod	Stoc k
CL03.0907.0500	500 ml	PE	0

CL03.2001

TISAB I Buffer - ISE buffer acc. NEN-ISO 15713:2011

NaCl 58 g + HAc 57 ml + NaCitr. 300 mg / l H₂O (pH 5.25 with NaOH)

Art. Nr.	Pack	Mod	Stoc k
CL03.2001.1000	1 l	PE	7
CL03.2001.5000	5 l	PE/HV	0
CL03.2001.9010	10 l	PE	0

CL12.2001

TISAB-ENOL Buffer - ISE buffer

10 g CDTA + 58 g NaCl + 57 ml HAc + 29.4 g NaCitr. / l H₂O to pH 5.5 with NaOH

Art. Nr.	Pack	Mod	Stoc k
CL12.2001.0250	250 ml	PE	0
CL12.2001.1000	1 l	PE	0

CL03.2002

TISAB II Buffer - ISE buffer

58 g NaCl + 57 ml CH₃COOH + 4 g CDTA / l H₂O (pH 5.3 with NaOH)

Art. Nr.	Pack	Mod	Stoc k
CL03.2002.0500	500 ml	PE	0
CL03.2002.1000	1 l	PE	3

CL03.2007

TISAB III Buffer - ISE buffer

265 g NH₄Cl + 53 g NH₄Ac + 19 g CDTA / l H₂O (pH 6 with NaOH) Do not store cold !

Art. Nr.	Pack	Mod	Stoc k
CL03.2007.0500	500 ml	PE	14
CL03.2007.1000	1 l	PE	3
	10 l	PE	0
CL03.2007.9010			

Art. Nr.	Pack	Mod	Stoc k
CL03.2009.0500	500 ml	PE	0

CL03.2010

TISAB TRIS Buffer - ISE buffer

121.14 g C₄H₁₁NO₃ (pH 9 with HNO₃)

Art. Nr.	Pack	Mod	Stock
----------	------	-----	-------

CL03.2010.1000	1 l	PE	0
----------------	-----	----	---

Art. Nr.	Pack	Mod	Stoc k
CL03.1416.0250	250 ml	PE	0

CL03.1415

Sodium chloride 58.4 g/l solution

58.4 g NaCl / l H₂O

Casnr.: 7647-14-5

Art. Nr.	Pack	Mod	Stoc k
CL03.1415.0250	250 ml	PE	0

CL03.1414

Sodium chloride 30 g/l solution

30 g NaCl / l H₂O

Casnr.: 7647-14-5

Art. Nr.	Pack	Mod	Stoc k
CL03.1414.0250	250 ml	PE	5

CL03.1413
Sodium chloride 23 g/l solution
 23 g NaCl / l H₂O
 Casnr.: 7647-14-5

Art. Nr.	Pack	Mod	Stoc k
CL03.1413.0250	250 ml	PE	0

CL03.1412
Sodium chloride 5.84 g/l solution
 5.84 g NaCl / l H₂O
 Casnr.: 7647-14-5

Art. Nr.	Pack	Mod	Stoc k
CL03.1412.0250	250 ml	PE	0

CL03.1411
Sodium chloride 3 g/l solution
 3 g NaCl / l H₂O
 Casnr.: 7647-14-5

Art. Nr.	Pack	Mod	Stoc k
CL03.1411.0250	250 ml	PE	0

По вопросам продаж и поддержки обращайтесь:

Алматы (727)345-47-04
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37
Пермь (342)205-81-47

Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Саранск (8342)22-96-24
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97
Тверь (4822)63-31-35

Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +(727)345-47-04

Беларусь +(375)257-127-884

Узбекистан +998(71)205-18-59

Киргизия +996(312)96-26-47

эл.почта: cmf@nt-rt.ru || сайт: <https://chem-lab.nt-rt.ru/>